






PATIENT INSTRUCTIONS

Kinex KA2™ Ankle CPM

To Change Extension/Flexion/Inversion/Eversion:

1. Press the **STOP** button.
2. Press the **MENU** button once; M1 should appear on the screen.
3. To change Extension, press the upper left **Extension** button on the hand controller. 
4. Press **(+)** button to increase and **(-)** button to decrease.
5. To change Flexion, press the inner left **Flexion** button. 
6. Press **(+)** button to increase and **(-)** button to decrease.
7. To change Inversion, press the inner right **Inversion** button. 
8. Press **(+)** button to increase and **(-)** button to decrease.
9. To change Eversion, press the upper right **Eversion** button. 
10. Press **(+)** button to increase and **(-)** button to decrease.
11. Press the **STOP** button to save the values.
12. Press the **START** button; device will go to midrange and stop.
13. Press the **START** button to begin therapy.



Distributed by:

KINEX
Phase I Rehab™