

PATIENT INSTRUCTIONS


Kinex KE2™ Elbow CPM

To change Extension or Flexion (straightening & bending of the elbow):

1. Press the **STOP** button.
2. Press the **STOP** button again and hold for 5 seconds.
3. An arrow will be pointing to the Extension value.
4. Press the plus (+) or minus (-) keys to go up or down.
Press the **Extension Flexion** button; arrow will point to Flexion.
5. Press the (+) or (-) keys to go up or down.

To change Pronation or Supination:

1. Press the **Pronation/Supination** button.
2. Arrow will be pointing to the Pronation value
3. Press the plus (+) or minus (-) keys to go up or down.
4. Press **Pronation/Supination** button.
5. Arrow will be pointing to the Supination value
6. Press the plus (+) or minus (-) keys to go up or down.
7. Press **START** again to recalibrate. CPM will go to midrange and stop.
8. Press **START** again.



Distributed by:

KINEX
Phase I Rehab™